


USER GUIDE



YELLOW BUTTON



TURN ON YOUR PINGCARES WATCH

- Hold for 3 seconds to power on
- The home screen will display time & date

RED BUTTON



EMERGENCY CALLING

- Press for 2 seconds to call for assistance
- Wait to connect or hit cancel if needed

6:08PM

Monday
Apr 29, 2024

PingCares

**SWIPE LEFT TO SEE
MORE SCREENS**



1289

Steps

**Check out the total
number of steps
you've walked today.**



Health



**"Tap" and
Measure
Vitals.**

Heart Rate

Blood Oxygen

Temperature



Contacts



**"Tap" and Find
Contacts to Call.**

MEASURING VITALS



**SWIPE TO HEALTH,
TAP ANYWHERE ON
THE SCREEN.**



**TAP ON THE VITAL YOU'D
LIKE TO MEASURE.**




**Make sure the watch
is worn securly on
your wrist in order to
measure your vitals.**



**HOLD STEADY AND
WAIT FOR THE WATCH
TO MEASURE.**

CALLING YOUR CONTACTS



SWIPE TO CONTACTS,
TAP ON THE SCREEN.



TAP "PAGE DOWN" TO
SEE MORE CONTACTS.



TAP ON YOUR
DESIRED CONTACT
TO CALL.



WAIT FOR THE CALL TO
CONNECT. TO HANG UP
ANYTIME, PRESS THE
YELLOW BUTTON

MENU SCREEN

SWIPE UP
TO GET TO
MENU SCREEN.



SWIPE DOWN
TO GET BACK TO
TIME SCREEN.



BATTERY
LEVEL



CELLULAR
STRENGTH



SCREEN
BRIGHTNESS



*Tap to adjust



LOWER
BRIGHTNESS



RINGER ON



*Tap to turn off



RINGER OFF

(This would SILENCE
all sounds)

USING THE CHARGER




To charge your watch, simply make sure it's plugged in and place it right on the stand! We recommend that you pick a set time, for example, breakfast time, everyday to charge your watch. Leave it on the charger

for about an hour, and it will be good to go for 24 hours.

Your watch will show the following statuses:



**Your watch
needs charging.**



**Charging is in
progress.**



**Your watch
is ready!**

WHAT DOES THIS MEAN?



Friendly reminder

In case your watch is forgotten somewhere, it will send you a friendly reminder to put it on!

How are you feeling?

Emotional health is important. If you see these faces on your watch, please tap on one that represents how you feel at the moment!



Happy



Neutral



Sad



FOUNDER STORY

Thank you for choosing PingCares. My name is Sean. I took this picture with my grandparents when I was a little boy. Despite them living in Beijing, I've always been close to my grandparents. However, I feel helpless at times when I can't be there should help be needed. During COVID, my worry was intensified and I felt a strong calling to provide peace of mind for family members who live far away from their elderly loved ones.



My background is in Engineering. I knew that technology could be a wonderful tool if it is made simple and practical. During grad school at the University of Chicago, I met many like-minded classmates. It became our mission to build a solution that improves the safety, wellness, and connection for our elders. Ping is named after my grandma, which is a Chinese word for peace and safety. Thank you for allowing us to serve you. We wish you great health and continued blessings.